



# EUROLINE-FOOD



## EUROLINE-FOOD Profile 216 (IgG) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
DP 3024-0404-1 G	Foods	IgG	Extract-coated membrane strips	4 x 04 (16)

**Indication:** The EUROLINE-FOOD test **supports** the diagnosis of sensitisations against foods and food additives, which may lead to unspecific health problems, by detection of IgG antibodies (e.g. gastrointestinal disorders, inflammatory skin diseases, migraine, chronic fatigue syndrome, etc.).

**Application:** The EUROLINE-FOOD test provides semiquantitative in vitro determination of the food-specific immunoglobulin class IgG in serum or plasma and contributes to the diagnosis of intolerances against foods and food additives. The test is a multi-parameter test which contains optimised combinations of relevant foods and their respective additives and therefore enables simultaneous analyses of IgG against both.

**Test principle:** The test kit contains test strips each coated with 54 different foods and food additives. The test strips are first moistened and then incubated with patient sample in the first reaction step. If samples contain specific antibodies of class IgG (and IgE, if present), they will bind to the antigenic components coated on the strip. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

### Contents of a test kit:

Description	Format	Symbol
1. <b>Test strip 1 coated with 54 foods and food additives and calibrators</b>	1 x 4 strips	
2. <b>Test strip 2 coated with 54 foods and food additives and calibrators</b>	1 x 4 strips	
3. <b>Test strip 3 coated with 54 foods and food additives and calibrators</b>	1 x 4 strips	
4. <b>Test strip 4 coated with 54 foods and food additives and calibrators</b>	1 x 4 strips	
5. <b>Enzyme conjugate</b> Alkaline phosphatase-labelled anti-human IgG (goat), ready for use	1 x 30 ml	
6. <b>Universal buffer</b> 10x concentrated	1 x 100 ml	
7. <b>Substrate solution</b> Nitroblue tetrazolium chloride/5-Bromo-4-chloro-3-indolylphosphate (NBT/BCIP), ready for use	1 x 30 ml	
8. <b>Incubation tray</b> , volume-reduced (400 µl)	2 x 10 channels	
9. <b>Plastic foil</b>	1 sheet	
10. <b>Instruction booklet</b>	1 booklet	
Lot description		Storage temperature
In vitro diagnostic medical device		Unopened usable until

Performance of the test requires incubation trays or other components, which are not provided in the test kits. They are available from EUROIMMUN under the following order numbers:

ZD 9897-0130 Incubation tray (volume-reduced 400 µl) with 30 channels (black, compatible with EUROBlotMaster and EUROBlotCamera System)

ZD 9897-0144 Incubation tray (volume-reduced 400 µl) with 44 channels (black, compatible with EUROBlotOne, EUROBlotMaster and EUROBlotCamera System)

**Modifications to the former version are marked in grey.**





# EUROLINE-FOOD



For the creation of work protocols and the evaluation of incubated test strips using **EUROLineScan** you require green paper and adhesive plastic foil:

ZD 9880-0101 Green paper (1 sheet)

ZD 9885-0116 Adhesive foil for approx. 16 test strips

ZD 9885-0130 Adhesive foil for approx. 30 test strips

For covering the incubation trays the adhesive foil can be used as well.

## Preparation and stability of the reagents

**Note:** All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. Unopened, reagents are stable until the indicated expiry date when stored at +2°C to +8°C. After initial opening, reagents are stable for 12 months or until the expiry date, unless stated otherwise below. Opened reagents must also be stored at +2°C to +8°C and protected from contamination.

- **Coated test strips:** Ready for use. Open the package with the test strips only when the strips have reached room temperature (+18°C to +25°C) to prevent condensation on the strips. After removal of the strips the package should be sealed tightly and stored at +2°C to +8°C.
- **Enzyme conjugate:** Ready for use. Mix thoroughly before using.
- **Universal buffer:** The universal buffer is supplied as a 10x concentrate. For the preparation of the working-strength universal buffer shake the bottle. The amount required should be removed from the bottle using a clean pipette tip and diluted 1:10 with deionised or distilled water. Due to the special membrane used for the present EUROLINE the working-strength universal buffer is used for the dilution of patient samples and the washing of the test strips. For the incubation of 1 test strip 5.0 ml buffer concentrate should be diluted with 45.0 ml water. The working-strength universal buffer should be used on the same working day.
- **Substrate solution:** Ready for use. Close bottle immediately after use, as the contents are sensitive to light ☀.

**Storage and stability:** The test kit must be stored at a temperature between +2°C and +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

**Waste disposal:** Patient samples and incubated test strips should be handled as infectious waste. Other reagents do not need to be collected separately, unless stated otherwise in official regulations.

**Warning:** Some of the reagents contain sodium azide in a non-declarable concentration. Avoid skin contact.

## Preparation and stability of the patient samples

**Sample material:** Human serum or EDTA, heparin or citrate plasma.

**Stability: Patient samples** to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

**Sample dilution:** The patient samples to be investigated are diluted 1:51 with the working-strength universal buffer using a clean pipette tip and mixed thoroughly by vortexing (for example, dilute 10 µl sample with 500 µl working-strength universal buffer). Sample pipettes are not suitable for mixing.



# EUROLINE-FOOD



## Incubation

**Pre-treatment:** For every sample place one test strip 1 to 4 in separate incubation trays (Make sure that the surface of the test strips is not damaged!). Fill each of the channels with 1.0 ml working-strength universal buffer and incubate the test strips for **5 minutes**. Afterwards aspirate off all the liquid.

**Sample incubation:** Manual:  
Fill each channel containing the test strip with **510 µl** of **1:51 diluted sample** using a clean pipette tip. Incubate for **60 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.

Automatic:  
For the automated version with **1:51 diluted sample** the incubation volume must be increased to **1020 µl** (1.0 ml working-strength universal buffer plus 20 µl sample).

**Washing:** Aspirate off the liquid from each channel and wash for **3 x 5 minutes** with 1.0 ml working-strength universal buffer on a rocking shaker.

**Conjugate incubation:** Pipette 1.0 ml enzyme conjugate (alkaline phosphatase-labelled anti-human IgG) into each channel and incubate for **60 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.

**Washing:** Aspirate off the liquid from each channel. Wash as described above.

**Substrate incubation:** Pipette 1.0 ml substrate solution into each channel and incubate for **10 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.

**Stopping:** Aspirate off the liquid from each channel and wash each test strip **3 x 1 minute** with deionised or distilled water.

**Evaluate:** Place the test strip on the evaluation protocol, air dry and evaluate.


For automated incubation with the **EUROBlotMaster** select the program **Euro 09 FOOD**.

For automated incubation with the **EUROBlotOne** select the program **EURO 09 FOOD**.



## EUROLINE-FOOD Profile 216 (IgG)

### Incubation protocol





# EUROLINE-FOOD



## Interpretation of results

**Handling:** After stopping the reaction using deionised or distilled water, place the incubated test strips onto the adhesive foil of the green work protocol (created beforehand in the EUROLineScan program) using a pair of tweezers. The position of the test strips can be corrected while they are wet. As soon as all test strips have been placed onto the protocol, they should be pressed hard using filter paper and left to air-dry. The drying process should take place without any direct light, in an environment as dark as possible. After they have dried, the test strips will be stuck to the adhesive foil. Incubated strips that are still moist show a background colouring that disappears when they are completely dry. Therefore the evaluation of the strips is only to take place after the strips have completely dried.

For **digital evaluation** follow the instructions in the EUROLineScan user manual. The code for entering the **test** into EUROLineScan is **FOOD strip 1** for test strip 1, **FOOD strip 2** for test strip 2, **FOOD strip 3** for test strip 3 and **FOOD strip 4** for test strip 4.


**Note:** EUROLineScan evaluates the test strips taking into account lot-specific specifications. These are transferred automatically to the software with the aid of barcodes on the EUROLINE. If the barcode is not readable, it is possible to enter it manually into EUROLineScan (see EUROLineScan instructions for use). The barcode number required can be found on the label of the folded card for the test strips.

**Caution:** A serum control band and 4 calibrators are located on the lower end of the strip. The correct performance of the incubation is confirmed after the automated evaluation of the staining intensity of the serum control band and the calibrators by means of the EUROLineScan program.

**Some samples might display a dark background staining of the membrane and white bands at the position of the antigens. These lighter bands should be interpreted as negative.**



# EUROLINE-FOOD




**Antigens:** The membrane strips are coated with the following foods and food additives.

EUROLINE-FOOD Profile 216 (IgG)				
Positio n	Strip 1 (3001-a)	Strip 2 (3001-b)	Strip 3 (3001-c)	Strip 4 (3001-d)
1	Barley flour	Cherry	Carob	Lime
3	Gluten	Grape, (white/ blue)	Rapeseed	Lychee
5	Oat flour	Kiwi	Duck meat	Mango
6	Rye flour	Lemon	Goat	Papaya
7	Spelt	Nectarine	Goose	Pomegranate
8	Wheat flour	Orange	Ostrich	Raspberry
9	Buckwheat flour	Pineapple	Quail	Rose hip
11	Flax seed	Strawberry	Rabbit	Cantaloupe
12	Corn	Watermelon	Roe deer	Raisin
13	Millet	Pear	Guinea fowl	Anise
15	Rice	Plum	Horse	Bay leaf
17	Beef	Grapefruit	Kefir	Chamomile
18	Chicken	Peach	Beta-Lactoglobulin	Caper
14	Lamb meat	Date fruit	Butter	Chive
15	Pork, cooked	Basil	Camembert	Clove
16	Turkey	Pepper, (black/ white)	Casein	Coriander
17	Cow's milk	Cinnamon	Emmental cheese	Cumin
18	Egg yolk	Garlic	Cottage cheese	Dill
19	Egg white	Mustard seed	Mozzarella	Ginger
20	Goat's cheese	Nutmeg	Processed cheese	Marjoram
21	Goat's milk	Oregano	Curd cheese	Saffron
22	Sheep's milk	Parsley	Bamboo shoots	Sage
23	Sheep's cheese	Peppermint	Brussel sprouts	Cayenne pepper
24	Yogurt	Poppy seed	Cauliflower	Curry
25	Eggplant	Rosemary	Chard	Tarragon
26	Beetroot	Thyme	Chinese cabbage	Hops
27	Bell pepper	Vanilla	Fennel	Mint
28	Broccoli	Almond	Gourd	Brazil nut
29	Carrot	Cashew nut	Jerusalem artichoke	Macadamia nut
30	Celery	Cocoa bean	Kale	Pine nut
31	Chili	Hazelnut	Radish	Chestnut
32	Cucumber	Peanut	Savoy cabbage	Cola nut
33	Horseradish	Pistachio	Sweet potato	Carp
34	Leek	Sesame	Vine leave	Squid
35	Olive	Sunflower seed	Green cabbage	Eel
36	Onion	Walnut	Shallot	Gilthead seabream
37	Potato	Coconut	Liquorice root	Haddock
38	Red cabbage	Mushroom mixture 1*	Snow pea	Pike
39	Tomato	Mushroom mixture 2*	Broad bean	Turbot
40	Turnip	Crayfish	Chickpea	Herring
41	Zucchini	Salmon	Mung bean	Lobster
42	Artichoke	Tuna	Kidney bean	Mackerel
43	Asparagus	Clam	Chicory	Octopus
44	Spinach	Prawn	Iceberg lettuce	Oyster
45	String bean	Anchovy	Rocket	Sardine
46	Pea	Swordfish	Avocado	Ocean perch
47	Soya bean	Trout	Blackberry	Sea bass
48	Lentil	Sole	Blueberry	Caviar
49	White bean	Codfish	Cranberry	Crab
50	Lettuce	Brewer's yeast	Red currant	Agar-agar
51	Corn salad	Baker's yeast	Black currant	Aloe vera
52	Apple	Honey	Fig	Green tea
53	Apricot	Coffee	Gooseberry	Baking powder
54	Banana	Black tea	Honeydew melon	Safflower oil
SC	Serum control band	Serum control band	Serum control band	Serum control band
C4	Calibrator 4	Calibrator 4	Calibrator 4	Calibrator 4
C3	Calibrator 3	Calibrator 3	Calibrator 3	Calibrator 3
C2	Calibrator 2	Calibrator 2	Calibrator 2	Calibrator 2
C1	Calibrator 1	Calibrator 1	Calibrator 1	Calibrator 1

\*Mushroom mixture 1: Oyster mushroom, White mushroom, Shiitake, Chanterelle  
Mushroom mixture 2: Bay Boletus, Boletus



# EUROLINE-FOOD



## Interpretation of results:

When using EUROLineScan the intensity of the bands is measured by the software and converted into classes 0 to 4. The specific IgG titers are also given in U/ml\*, which are calculated for each test strip in relation to the intensity of the calibrators.

The classes can be divided into the following concentrations:

Class	Concentration [U/ml*]	Result	Interpretation
0	$0 < \text{IgG} \leq 7.5$	Negative	No specific antibodies detectable
1	$7.5 < \text{IgG} \leq 12.5$	Weak	Antibody titer in physiological range, no disease association
2	$12.5 < \text{IgG} \leq 20$	Increased	Sensitisation possible
3	$20 < \text{IgG} \leq 50$	Positive	Indicating sensitisation, food intolerance likely
4	$> 50$	Strongly positive	Very high antibody titer, almost always associated with sensitisation and manifest intolerance

\*The EUROLINE-FOOD tests were calibrated using the WHO reference serum “1<sup>st</sup> IRP 67/86”. The measurement values are given in units per milliliter [U/ml]. The limits given are reference values obtained in internal analysis.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

## Test characteristics

**Measurement range:** The EUROLINE is a semiquantitative method. The measurement range is given in classes 0 to 4.

**Limitations of the determination:** The antibody determination indicates the patient's level of sensitisation to the different foods and food additives investigated. It does not, however, provide information on the connection between the height of the IgG antibody titer determined with the EUROLINE-FOOD and the occurrence or severity of food intolerance.

Furthermore, it cannot be excluded that epitopes become inaccessible for binding to the antibodies during coating of the strips, which can result in missing signals and thus in false negatives. It is also possible that antibody binding sites are only created or become free for binding during (industrial) processing or the digestion process. These cases can also cause false negative results.

**Cross-reactions:** If the antigenic structures are similar, e.g. in chemically similar substances or due to a close botanic relationship, cross-reactions between IgG antibodies of different foods and food additives may occur. The specific IgG antibodies that have developed in a patient also attach to identical epitopes of a homologous antigen.

**Interferences:** Haemolytic, lipaemic and icteric samples up to a concentration of 5 mg/ml for haemoglobin, of 20 mg/ml for triglycerides and of 0.4 mg/ml for bilirubin showed no effect on the analytical results of the EUROLINE-FOOD.

**Inter- and intra-assay variation:** The inter-assay variation was determined by multiple analyses of characterised samples over several days. The intra-assay variation was determined by multiple analyses of characteristic samples on one day. In every case, the intensity of the bands was within the specified range. The EUROLINE-FOOD displays an excellent inter-assay and intra-assay reproducibility.





## Clinical significance

Food intolerance (food hypersensitivity) is an immunological reaction to food or food additives. The immune response arises due to an increased permeability of the small intestine, which allows food components to leak through the intestinal wall and into the bloodstream. This permeability (so-called leaky gut syndrome) can be caused by diet-related hyperacidity of the gut flora, medication, infections, preservatives, alcohol/nicotine, stress and further factors.

The immune system identifies these antigenic structures as foreign and starts to produce specific antibodies of class IgG (sensitisation). This immune response, which typically occurs after a certain period of time (several hours to days after intake of the food), triggers an inflammatory reaction, which may become chronic. Chronic inflammation can manifest in the entire body, not only in the gastrointestinal tract. The most frequent symptoms are diarrhoea, constipation, bloating, nausea, upset stomach, irritable colon, migraine, asthma, joint disorders, lack of concentration, skin disorders and weight problems (over-/underweight). To relieve symptoms, an elimination diet (exclusion diet) for a set period of time is usually recommended. Those foods for which (high) IgG antibody concentrations were measured are excluded from the patient's diet. Improving the gut flora with therapeutic measures can also be taken into account to reduce or prevent the permeability of the intestinal wall to food antigens.

The association with different diseases has been described in many studies. It is often accompanied by high levels of IgG antibodies against various foods or food additives. This is the case in patients with headaches/migraine, arthritis, irritable colon (irritable bowel syndrome), inflammatory bowel diseases, asthma and a number of other diseases.

**Headaches and migraine:** A study from the 1930s already showed that 66% of 127 migraine patients were free of complaints after undertaking an elimination diet (Sheldon and Randolph, 1935). In other more recent studies these results were confirmed and it was shown that an elimination diet based on IgG levels significantly influences the frequency of migraine attacks (Alpay et al., 2010; Arroyave Hernández et al., 2007; Aydinlar et al., 2012; Mitchell et al., 2011).

**Arthritis:** Food intolerance has been known for more than 35 years to be a causative factor for developing arthritis. In a study with 22 patients with rheumatoid arthritis undertaking an elimination diet, 20 patients (91%) experienced an improvement of symptoms. 19 patients again suffered a deterioration of their health status when the reactive foods were reintroduced into their diet (Hicklin et al., 1980).

**Irritable colon:** In various studies with patients with irritable bowel syndrome the participants were treated with an elimination diet. A significant number of patients experienced an improvement of symptoms up to an at least subjective full recovery (Jones et al., 1982; Nanda et al., 1989). Atkinson and colleagues, and Drisko and colleagues performed studies with respect to IgG-based elimination diet in irritable-colon patients. Both groups showed that the symptoms could be reduced significantly (Atkinson et al., 2004; Drisko et al., 2006).

**Intestinal bowel diseases:** In a retrospective study Cai and colleagues (2014) investigated 112 patients with inflammatory bowel diseases [Crohn's disease (n = 79) and ulcerative colitis (n = 33)] and compared the data to those obtained in a healthy control group (n = 266). They found out that patients with inflammatory bowel diseases have a high prevalence of antibodies of class IgG against certain foods and that this can be used for establishing an elimination diet (Cai et al., 2014). Further studies yielded similar results (Lindberg et al., 1992; Bentz et al., 2010; Kawaguchi et al., 2014).

Increased levels of IgG antibodies against foods or food additives have also been observed in **asthma** (Calderon et al., 2010; Codina et al., 1997), **autism** (Lucarelli et al., 1995; Kidd, 2002) or **atopic dermatitis** (Shakib et al., 1977).

The internet-based study of Gaby (1998) gives an overview of the possible association of food intolerance with various diseases.





The diagnostic value of IgG (or IgG4) antibody detection for the diagnosis of food intolerance is being controversially discussed. Different studies with up to several thousand patient samples, however, have come to the conclusion that the determination of the antibody titer can be a useful tool for the identification of food intolerances and for targeted patient therapy (Atkinson et al., 2004; Zar et al., 2005; Bernardi et al., 2008; Volpi and Maccari 2009). Generally, it could be shown that elimination of those foods, against which strong antibody reactions were measured, helped to improve symptoms or to promote complete recovery in a statistically significant number of patients.

In a large internet-based study, Mullin and colleagues investigated the diagnostic benefit of IgG determination as a marker of food intolerance, and its clinical relevance. The authors concluded that IgG determination provides a clinically useful basis for establishing an elimination diet (Mullin et al., 2010).

## Literature references

1. Alpay K, Ertas M, Orhan EK, Ustay DK, Lieners C, Baykan B. **Diet restriction in migraine, based on IgG against foods: a clinical double-blind, randomised, cross-over trial.** Cephalalgia. 2010 Jul;30(7):829-37. doi: 10.1177/0333102410361404. Epub 2010 Mar 10.
2. Arroyave Hernández CM, Echavarría Pinto M, Hernández Montiel HL. **Food allergy mediated by IgG antibodies associated with migraine in adults.** Rev Alerg Mex. 2007 Sep-Oct;54(5):162-8. Erratum in: Rev Alerg Mex. 2010 Mar-Apr;57(2):49. Echevarría Pinto, Mauro.
3. Awazuhara H, Kawai H, Maruchi N. **Major allergens in soybean and clinical significance of IgG4 antibodies investigated by IgE- and IgG4-immunoblotting with sera from soybean-sensitive patients.** Clin Exp Allergy. 1997 Mar;27(3):325-32.
4. Aydinlar EI, Dikmen PY, Tiftikci A, Saruc M, Aksu M, Gunsoy HG, Tozun N. **IgG-based elimination diet in migraine plus irritable bowel syndrome.** Headache. 2013 Mar;53(3):514-25. doi: 10.1111/j.1526-4610.2012.02296.x. Epub 2012 Dec 6.
5. Bentz S, Hausmann M, Piberger H, Kellermeier S, Paul S, Held L, Falk W, Obermeier F, Fried M, Schölmerich J, Rogler G. Digestion. **Clinical relevance of IgG antibodies against food antigens in Crohn's disease: a double-blind cross-over diet intervention study.** 2010;81(4):252-64. doi: 10.1159/000264649. Epub 2010 Jan 30.
6. Bernardi D, Borghesan F, Faggian D, Bianchi FC, Favero E, Billeri L, Plebani M. **Time to reconsider the clinical value of immunoglobulin G4 to foods?** Clin Chem Lab Med. 2008;46(5):687-90.
7. Bischoff SC, Herrmann A, Manns MP. **Prevalence of adverse reactions to food in patients with gastrointestinal disease.** Allergy. 1996 Nov;51(11):811-8.
8. Cai C, Shen J, Zhao D, Qiao Y, Xu A, Jin S, Ran Z, Zheng Q. **Serological investigation of food specific immunoglobulin G antibodies in patients with inflammatory bowel diseases.** PLoS One. 2014 Nov 13;9(11):e112154. doi: 10.1371/journal.pone.0112154. eCollection 2014.
9. Calderon TE, Ferrero M, Marino GM, Cordoba A, Beltramo D, Muino JC, Rabinovich GA, Romero MD. **Meat-specific IgG and IgA antibodies coexist with IgE antibodies in sera from allergic patients: clinical association and modulation by exclusion diet.** J Biol Regul Homeost Agents. 2010 Jul-Sep;24(3):261-71.
10. Codina RM, Calderón E, Lockey RF, Fernández-Caldas E, Rama R. **Specific immunoglobulins to soybean hull allergens in soybean asthma.** Chest. 1997 Jan;111(1):75-80.
11. Drisko J, Bischoff B, Hall M, McCallum R. **Treating irritable bowel syndrome with a food elimination diet followed by food challenge and probiotics.** J Am Coll Nutr. 2006 Dec;25(6):514-22.
12. Duchateau J, Michils A, Lambert J, Gossart B, Casimir G. **Anti-betalactoglobulin IgG antibodies bind to a specific profile of epitopes when patients are allergic to cow's milk proteins.** Clin Exp Allergy. 1998 Jul;28(7):824-33.



13. el Rafei A, Peters SM, Harris N, Bellanti JA. **Diagnostic value of IgG4 measurements in patients with food allergy.** Ann Allergy. 1989 Feb;62(2):94-9.
14. Gaby AR. **The role of hidden food allergy/intolerance in chronic disease.** Altern Med Rev. 1998 Apr;3(2):90-100.
15. Germano P, Pezzini A, Boccagni P, Zanoni G, Tridente G. **Specific humoral response to cows' milk proteins and ovalbumin in children with atopic dermatitis.** Int J Clin Lab Res. 1993;23(4):206-11.
16. Haddad ZH, Vetter M, Friedmann J, Sainz C, Brunner E. **Detection and kinetics of antigen-specific IgE and IgG immune complexes in food allergy.** Ann Allergy. 1983 Aug;51(2 Pt 2):255.
17. Halpern GM, Scott JR. **Non-IgE antibody mediated mechanisms in food allergy.** Ann Allergy. 1987 Jan;58(1):14-27. Review.
18. Jones VA, McLaughlan P, Shorthouse M, Workman E, Hunter JO. **Food intolerance: a major factor in the pathogenesis of irritable bowel syndrome.** Lancet. 1982 Nov 20;2(8308):1115-7.
19. Kawaguchi T, Mori M, Saito K, Suga Y, Hashimoto M, Sako M, Yoshimura N, Uo M, Danjo K, Ikenoue Y, Oomura K, Shinozaki J, Mitsui A, Kajiura T, Suzuki M, Takazoe M. **Food antigen-induced immune responses in Crohn's disease patients and experimental colitis mice.** J Gastroenterol. 2015 Apr;50(4):394-405. doi: 10.1007/s00535-014-0981-8. Epub 2014 Aug 7.
20. Kidd PM. **Autism, an extreme challenge to integrative medicine.** Part 2: medical management. Altern Med Rev. 2002 Dec;7(6):472-99.
21. Lindberg E, Magnusson KE, Tysk C, Järnerot G. **Antibody (IgG, IgA, and IgM) to baker's yeast (*Saccharomyces cerevisiae*), yeast mannan, gliadin, ovalbumin and betalactoglobulin in monozygotic twins with inflammatory bowel disease.** Gut. 1992 Jul;33(7):909-13.
22. Lucarelli S, Frediani T, Zingoni AM, Ferruzzi F, Giardini O, Quintieri F, Barbato M, D'Eufemia P, Cardi E. **Food allergy and infantile autism.** Panminerva Med. 1995 Sep;37(3):137-41.
23. Mitchell N, Hewitt CE, Jayakody S, Islam M, Adamson J, Watt I, Torgerson DJ. **Randomised controlled trial of food elimination diet based on IgG antibodies for the prevention of migraine like headaches.** Nutr J. 2011 Aug 11;10:85. doi: 10.1186/1475-2891-10-85.
24. Mullin GE, Swift KM, Lipski L, Turnbull LK, Rampertab SD. **Testing for food reactions: the good, the bad, and the ugly.** Nutr Clin Pract. 2010 Apr;25(2):192-8. doi: 10.1177/0884533610362696.
25. Nanda R, James R, Smith H, Dudley CR, Jewell DP. **Food intolerance and the irritable bowel syndrome.** Gut. 1989 Aug;30(8):1099-104.
26. Shakib F, Brown HM, Phelps A, Redhead R. **Study of IgG sub-class antibodies in patients with milk intolerance.** Clin Allergy. 1986 Sep;16(5):451-8.
27. Shakib F, McLaughlan P, Stanworth DR, Smith E, Fairburn E. **Elevated serum IgE and IgG4 in patients with atopic dermatitis.** Br J Dermatol. 1977 Jul;97(1):59-63.
28. Sheldon JM, Randolph TG. **Allergy in migraine-like headaches.** Am J Med Sci 1935;190:232-236.
29. Volpi N, Maccari F. **Serum IgG responses to food antigens in the italian population evaluated by highly sensitive and specific ELISA test.** J Immunoassay Immunochem. 2009;30(1):51-69. doi: 10.1080/15321810802571903.
30. Zar S, Benson MJ, Kumar D. **Food-specific serum IgG4 and IgE titers to common food antigens in irritable bowel syndrome.** Am J Gastroenterol. 2005 Jul;100(7):1550-7.



---

# EUROLINE-FOOD

---





---

# EUROLINE-FOOD

---



DP\_3024-1\_A\_UK\_C04.doc  
Version: 07/12/2018

