

















## Anti-SS-A ELISA (IgG) Test instruction

| ORDER NO.      | ANTIBODIES AGAINST | IG CLASS | SUBSTRATE                  | FORMAT       |
|----------------|--------------------|----------|----------------------------|--------------|
| EA 1595-9601 G | SS-A               | IgG      | Ag-coated microplate wells | 96 x 01 (96) |

**Indications:** Sjögren's syndrome, systemic lupus erythematosus.

**Principles of the test:** The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human autoantibodies of the IgG class against SS-A in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with SS-A. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

### Contents of the test kit:

| Component   | Colour  | Format  | Symbol  |
|---|---|---|---|
| <b>1. Microplate wells</b> ,<br>coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use | ---   | 12 x 8  |                          |
| <b>2. Calibrator 1</b><br>200 RU/ml (IgG, human), ready for use   | dark red  | 1 x 2.0 ml  |                         |
| <b>3. Calibrator 2</b><br>20 RU/ml (IgG, human), ready for use  | red   | 1 x 2.0 ml  |                        |
| <b>4. Calibrator 3</b><br>2 RU/ml (IgG, human), ready for use   | light red   | 1 x 2.0 ml  |                        |
| <b>5. Positive control</b><br>(IgG, human), ready for use   | blue  | 1 x 2.0 ml  |                        |
| <b>6. Negative control</b><br>(IgG, human), ready for use   | green   | 1 x 2.0 ml  |                        |
| <b>7. Enzyme conjugate</b><br>peroxidase-labelled anti-human IgG (rabbit), ready for use  | green   | 1 x 12 ml   |                        |
| <b>8. Sample buffer</b><br>ready for use  | light blue  | 1 x 100 ml  |                        |
| <b>9. Wash buffer</b><br>10x concentrate  | colourless  | 1 x 100 ml  |                        |
| <b>10. Chromogen/substrate solution</b><br>TMB/H <sub>2</sub> O <sub>2</sub> , ready for use  | colourless  | 1 x 12 ml   |                        |
| <b>11. Stop solution</b><br>0.5 M sulphuric acid, ready for use   | colourless  | 1 x 12 ml   |                        |
| <b>12. Test instruction</b>   | ---   | 1 booklet   |   |
| <b>13. Quality control certificate</b>  | ---   | 1 protocol  |   |
|  Lot   |  |  Storage temperature |  Unopened usable until |
|  In vitro determination  |   |   |   |


**Storage and stability:** The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

**Waste disposal:** Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.



## Preparation and stability of the reagents

**Note:** All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).  
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).  
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.  
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light . The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

**Warning:** Calibration and controls used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays or indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

## Preparation and stability of the patient samples

**Sample material:** Human serum or EDTA, heparin or citrate plasma.

**Stability: Patient samples** to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

**Sample dilution: Patient samples** are diluted **1:201** in sample buffer. For example: dilute 5 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.



## Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

### (Partly) manual test performance

**Sample incubation:** (1<sup>st</sup> step) Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

**Washing:** Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.  
Automatic: Wash reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

**Conjugate incubation:** (2<sup>nd</sup> step) Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

**Washing:** Empty the wells. Wash as described above.

**Substrate incubation:** (3<sup>rd</sup> step) Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) protect from direct sunlight.

**Stopping the reaction:** Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

**Measurement:** **Photometric measurement** of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



## Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry.

Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.

## Pipetting protocol

|   | 1    | 2    | 3    | 4    | 5 | 6 | 7    | 8    | 9    | 10   | 11 | 12 |
|---|------|------|------|------|---|---|------|------|------|------|----|----|
| A | C 2  | P 6  | P 14 | P 22 |   |   | C 1  | P 4  | P 12 | P 20 |    |    |
| B | pos. | P 7  | P 15 | P 23 |   |   | C 2  | P 5  | P 13 | P 21 |    |    |
| C | neg. | P 8  | P 16 | P 24 |   |   | C 3  | P 6  | P 14 | P 22 |    |    |
| D | P 1  | P 9  | P 17 |      |   |   | pos. | P 7  | P 15 | P 23 |    |    |
| E | P 2  | P 10 | P 18 |      |   |   | neg. | P 8  | P 16 | P 24 |    |    |
| F | P 3  | P 11 | P 19 |      |   |   | P 1  | P 9  | P 17 |      |    |    |
| G | P 4  | P 12 | P 20 |      |   |   | P 2  | P 10 | P 18 |      |    |    |
| H | P 5  | P 13 | P 21 |      |   |   | P 3  | P 11 | P 19 |      |    |    |

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. Therefore, the number of tests performed can be matched to the number of samples, minimizing reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

## Calculation of results


**Semiquantitative:** Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

|                       |                 |
|-----------------------|-----------------|
| <b>Ratio &lt;1.0:</b> | <b>negative</b> |
| <b>Ratio ≥1.0:</b>    | <b>positive</b> |

**Quantitative:** The standard curve from which the concentration of antibodies in the serum samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction of a serum sample lies above the value of calibrator 1 (200 RU/ml). The result should be given as “>200 RU/ml“. It is recommended that the sample be re-tested at a dilution of 1:800. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range (**cut-off**) recommended by EUROIMMUN is 20 relative units (RU)/ml. EUROIMMUN recommends interpreting results as follows:

|                      |                 |
|----------------------|-----------------|
| <b>&lt;20 RU/ml:</b> | <b>negative</b> |
| <b>≥20 RU/ml:</b>    | <b>positive</b> |

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

## Test characteristics

**Calibration:** As no international reference serum exists for the quantitative detection of antibodies against SS-A, the calibration is performed in relative units (RU)/ml. The reactivity of the Anti-SS-A ELISA was verified using the human reference serum CDC-ANA #7 of the “Centers for Disease Control“ (Atlanta, USA).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result. minutes



**Antigen:** The microplate wells were coated with SS-A purified by affinity chromatography from calf thymus.

The SS-A antigen is localized in the cell nucleus and is involved in the processing of mRNA to translationally active molecules. It is a small ribonucleoprotein consisting of an RNA molecule (Y1, Y2, Y3, Y4 or Y5 RNA: 80-112 bases in length) and two different proteins, which are the targets for antibodies against SS-A. Initially, in 1984, a protein of 60 kDa was described as a component of ribonucleoproteins. In 1988, antibodies against a further protein with a molecular weight of 52 kDa (Ro-52) were detected in anti-SS-A-positive sera by means of Westernblot. However, Ro-52 does not appear to be a stable component of the native ribonucleoprotein particle.

Only test systems which use native SS-A 60 kDa as the antigen should be used for the detection of autoantibodies against SS-A in SLE or Sjögren's syndrome (SS). Test systems which include Ro-52 in the antigen substrate are not recommended, since antibodies against Ro-52 are frequently also found in myositis patients. Thus, the specificity of these test systems for SLE and SS is reduced. All anti-SS-A positive sera from SLE or SS patients can be identified using native SS-A 60 kDa as the antigen.

**Linearity:** The linearity of the Anti-SS-A ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination  $R^2$  for all sera was  $> 0.95$ . The Anti-SS-A ELISA (IgG) is linear at least in the tested concentration range (13 RU/ml to 185 RU/ml).

**Detection limit:** : The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-SS-A ELISA (IgG) is 1 RU/ml.

**Cross reactivity:** This ELISA showed no cross reactivity.

**Interference:** Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for hemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

**Reproducibility:** The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

| <i>Intra-assay variation, n = 20</i> |                    |        |
|--------------------------------------|--------------------|--------|
| Serum                                | Mean value (RU/ml) | CV (%) |
| 1                                    | 92                 | 3.0    |
| 2                                    | 113                | 1.9    |
| 3                                    | 114                | 4.2    |

| <i>Inter-assay variation, n = 4 x 6</i> |                    |        |
|---|--------------------|--------|
| Serum                                   | Mean value (RU/ml) | CV (%) |
| 1                                       | 92                 | 3.1    |
| 2                                       | 116                | 4.1    |
| 3                                       | 121                | 3.0    |

**Reference range:** The levels of the anti-SS-A antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 206 healthy blood donors. With a cut-off of 20 RU/ml, all blood donors were anti-SS-A negative.



## Clinical significance

Antibodies (AAb) against nuclear antigens (ANA) are directed against various cell nuclear components. Among the most important nuclear antigens, including cytoplasmic antigens, are nRNP/Sm, Sm, SS-A (Ro), SS-B (La), Scl-70, PM-Scl, Jo-1, centromeres, PCNA, dsDNA, nucleosomes, histones and ribosomal P-proteins. They are mainly components of functional nuclear particles, are bound to nucleic acids or fulfil functions in the cell cycle, e.g. in transcription or translation.

The investigation of ANA and subsequent differentiation within the ANA (or ENA) spectrum contributes greatly to establishing a diagnosis, particularly in the following rheumatic diseases:

- systemic lupus erythematosus (SLE),
- Sharp syndrome (mixed connective tissue disease = MCTD),
- Sjögren's syndrome (SS),
- systemic sclerosis (SSc), and
- poly-/dermatomyositis (PM/DM).

**Sjögren's syndrome (SS)** is a chronic inflammatory autoimmune disease of the exocrine glands which can be found in one to four million people in the US alone. Nine out of ten patients are women. The main clinical feature of primary SS is ocular and oral dryness as a result of the destruction of lachrymal and salivary glands by lymphocytic infiltration. The pancreatic glands, the mucous secreting glands of the intestine, bronchia or vagina and the sudoriferous glands may also be affected. Around 5% of SS patients develop malignant lymphoma. In secondary SS, primary SS symptoms accompany rheumatoid arthritis (RA), SSc, SLE, PM/DM, primary biliary cirrhosis or autoimmune hepatitis.

**Anti-SS-A** are detected in 40%-95% of SS cases. They mostly occur in parallel with autoantibodies against SS-B (anti-La). Autoantibodies against SS-A are also found in 20%-60% of SLE patients and in neonatal lupus erythematosus (neonatal LE syndrome) with a prevalence of 100%. The antibodies are transmitted diaplacentally to the foetus and often cause congenital AV block in addition to inflammatory reactions when the mother is anti-SS-A or anti-SS-B positive (level I-III).

Note: Differentiation of anti-SS-A antibodies from those against the so-called Ro52 antigen (52 kDa protein, RING dependent E3 ligase) is of decisive diagnostic importance, since antibodies against Ro52 are not disease-specific, but are also detected in myositis, systemic sclerosis, neonatal lupus erythematosus and other collagenoses, primary biliary cirrhosis, autoimmune hepatitis and viral hepatitis.

| Antibodies against | Disease                            | Prevalence |
|--------------------|------------------------------------|------------|
| SS-A (Ro)          | Sjögren's syndrome (SS)            | 40% - 95%  |
|                    | Systemic lupus erythematosus (SLE) | 20% - 60%  |
|                    | Neonatal lupus erythematosus       | 95% - 100% |

## Literature references

1. Alba P, Bento L, Cuadrado MJ, Karim Y, Tungekar MF, Abbs I, Khamashta MA, D'Cruz D, Hughes GR. **Anti-dsDNA, anti-Sm antibodies, and the lupus anticoagulant: significant factors associated with lupus nephritis.** Ann Rheum Dis 62 (2003) 556-560.
2. Boire G, Gendron M, Monast N, Bastin B, Ménard HA. **Purification of antigenically intact Ro ribonucleoproteins; biochemical and immunological evidence that the 52-kD protein is not a Ro protein.** Clin Exp Immunol 100 (1995) 489-498.
3. EUROIMMUN AG. Stöcker W, Schlumberger W, Krüger C. **Alle Beiträge zum Thema Autoimmundiagnostik.** In: Gressner A, Arndt T (Hrsg.) Lexikon der Medizinischen Laboratoriumsdiagnostik. 2. Auflage. Springer Medizin Verlag, Heidelberg (2012).
4. Hartung K, Seelig HP. **Laboratory diagnostics of systemic autoimmune diseases. Part 1. Collagenoses.** [Article in German] Z Rheumatol 65 (2006) 709-724.
5. Ho KT, Reveille JD. **The clinical relevance of autoantibodies in scleroderma.** Arthritis Res Ther 5 (2003) 80-93.
6. Meyer\* W, Scheper\* T, Wilhelm\* K, Jarzabek-Chorzelska M, Kolacinska-Strasz Z, Schlumberger\* W, Stöcker\* W. (\*EUROIMMUN AG) **Antibodies against SS-A can only be precisely detected using the native antigen: Results of a study using the EUROLINE-WB.** In: Conrad K et al. (Hrsg). From Proteomics to Molecular Epidemiology: Relevance of Autoantibodies. Pabst Science Publishers 3: (2002) 651-652.
7. Phan TG, Wong RC, Adelstein S. **Autoantibodies to extractable nuclear antigens: making detection and interpretation more meaningful.** Clin Diagn Lab Immunol 9 (2002) 1-7.
8. Tan EM, Chan EK, Sullivan KF, Rubin RL. **Antinuclear antibodies (ANAs): diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity.** Clin Immunol Immunopathol 47 (1988) 121-141.
9. Tzioufas AG, Voulgarelis M. **Update on Sjögren's syndrome autoimmune epithelitis: from classification to increased neoplasias.** Best Pract Res Clin Rheumatol 21 (2007) 989-1010.
10. Vitali C, Bombardieri S, Jonsson R, Moutsopoulos HM, Alexander EL, Carsons SE, Daniels TE, Fox PC, Fox RI, Kassan SS, Pillemer SR, Talal N, Weisman MH, and the European Study Group on Classification Criteria for Sjögren's Syndrome. **Classification criteria for Sjögren's syndrome: a revised version of the European criteria proposed by the American-European Consensus Group.** Ann Rheum Dis 61 (2002) 554-558.