

Anti-GBM ELISA (IgG) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EA 1251-9601 G	GBM	IgG	Ag-coated microplate wells	96 x 01 (96)

Indication: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human autoantibodies of the immunoglobulin class IgG against antigens of the glomerular basement membrane (GBM antigens) in serum or plasma for the diagnosis of anti-GBM disease (Goodpasture's syndrome).

Application: The detection of circulating autoantibodies against glomerular basement membrane (GBM) is an important contribution to the diagnosis of anti-GBM disease, a progressive glomerulonephritis that often involves the lungs (then called Goodpasture syndrome). In addition to securing a diagnosis, the anti-GBM ELISA (IgG) is used for therapy monitoring. High autoantibody titers indicate an unfavourable development.

Principles of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with GBM antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.


Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	STRIPS
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	CAL 1
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	CAL 2
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	CAL 3
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	POS CONTROL
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	NEG CONTROL
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	CONJUGATE
8. Sample buffer ready for use	light blue	1 x 100 ml	SAMPLE BUFFER
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	WASH BUFFER 10x
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	SUBSTRATE
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	STOP SOLUTION
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
<div style="display: flex; justify-content: space-between; align-items: center;"> <div style="display: flex; gap: 10px;"> <div style="border: 1px solid black; padding: 2px;">LOT</div> Lot description</div> <div style="font-size: 2em; font-weight: bold;">CE</div> <div style="display: flex; gap: 10px;"> <div style="font-size: 1.2em;">⚠</div> Storage temperature</div> </div>			
<div style="display: flex; justify-content: space-between; align-items: center;"> <div style="display: flex; gap: 10px;"> <div style="border: 1px solid black; padding: 2px;">IVD</div> In vitro diagnostic medical device</div> <div style="font-size: 1.2em;">🕒</div> Unopened usable until</div>			



Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to +37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light . The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer. For example: dilute 10 µl sample in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The calibrators and controls are prediluted and ready for use, do not dilute them.



Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step) Transfer 100 µl of the calibrators, positive or negative control or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Mode").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (>10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step) Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Empty the wells. Wash as described above.

Substrate incubation: (3rd step) Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) (protect from direct sunlight).

Stopping: Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement: **Photometric measurement** of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.



Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open-system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1 to 4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7 to 10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimises reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <1.0:	negative
Ratio ≥1.0:	positive



Quantitative: The standard curve from which the concentration of antibodies in the serum samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as ">200 RU/ml". It is recommended that the sample be retested at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range (**cut-off**) recommended by EUROIMMUN is 20 relative units (RU)/ml. EUROIMMUN recommends interpreting results as follows:

<20 RU/ml:	negative
≥20 RU/ml:	positive

For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends retesting the samples.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against GBM antigens, the calibration is performed in relative units (RU/ml).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio values determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The microplate wells were coated with highly purified alpha-3-chains of collagen type IV which contain the relevant antigenic epitopes of the GBM antigen.



Main components of the glomerular basement membrane (GBM) are the extracellular matrix proteins collagen type IV, laminin, fibronectin and proteoglycans. Autoantibodies against the glomerular basement membrane target epitopes on the collagen type IV. The collagen type IV molecule consists of 3 alpha-chains with molecular weights of 170 kDa each. The alpha-chains form several triple helical domains separated by sequences incompatible with helix formation. A compact helical zone (7S domain) is located at the amino-terminal end and a globular knob (NC1 domain) at the carboxy-terminal end. The major antigen of anti-GBM antibodies is the NC1 domain of the alpha-3(IV)-chain.

Linearity: The linearity of the Anti-GBM ELISA (IgG) was determined by assaying at least 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was >0.95 . The Anti-GBM ELISA (IgG) is linear at least in the tested concentration range (15 RU/ml to 200 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-GBM ELISA (IgG) is 1.6 RU/ml.

Cross reactivity: This ELISA specifically detects autoantibodies of class IgG against GBM. Cross reactions with other autoantibodies were not found in samples from patients with the following antibodies: M2, MPO and PR3 (n = 15).

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 samples. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed on 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Sample	Mean value (RU/ml)	CV (%)
1	58	3.3
2	60	4.9
3	83	3.4

<i>Inter-assay variation, n = 4 x 6</i>		
Sample	Mean value (RU/ml)	CV (%)
1	60	5.1
2	65	6.6
3	85	5.0

Sensitivity and specificity: Sera from 14 patients with Goodpasture's syndrome and 180 patients with other diseases were examined with the EUROIMMUN Anti-GBM ELISA. The test showed a sensitivity and a specificity of 100% each.

Patient group (n = 194)	n	anti-GBM positive
Goodpasture's syndrome	14	14 (100%)
SLE	100	0
Wegener's granulomatosis	50	0
Microscopic polyangitis	20	0
Churg-Strauss syndrome	10	0

Reference range: The levels of the anti-GBM antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 350 healthy blood donors. With a cut-off of 20 RU/ml, all blood donors were anti-GBM negative.



Correlation of the ELISA with the indirect immunofluorescence test (IIFT): 14 sera from patients with Goodpasture's syndrome and 420 uncharacterized sera were investigated with the EUROIMMUN Anti-GBM ELISA (IgG) and the EUROIMMUN IIFT (IgG) was used as a reference. The test showed a specificity of 99% and a sensitivity of 100% with regard to the IIFT.

n = 434		IIFT	
		positive	negative
ELISA	positive	47	5
	negative	0	382

Clinical significance

Autoantibodies against glomerular basement membrane (GBM) can be serologically detected in Anti-GBM disease, a glomerular nephritis which frequently involves the lungs (then called Goodpasture syndrome). The prevalence of GBM antibodies in patients with lung involvement amounts to 80-90%, in patients without lung involvement to 60%. Approximately 10% of the patients with GBM antibodies also show additional MPO ANCA. The epitope of GBM antibodies is the non-collagen domain 1 (NC 1 domain) of collagen type IV. The collagen is expressed at the glomerular basement membrane of the kidney and the alveolar basal membrane of the lungs. The autoantibodies are directed against the glomerular as well as the alveolar basal membranes, which explains the specific organ manifestations of the disease. The autoantibodies activate the complement cascade with the release of chemoactive fragments and infiltration of the tissue through inflammatory cells such as monocytes/macrophages and neutrophils. In addition to the antibody production, also the T-cell mediated immune response plays a decisive role in the pathogenesis. While T helper cells help to initiate the antibody production, cytotoxic T cells are involved in the organ damage of kidneys and lungs. A genetic association of the disease exists to the alleles HLA-DRB1(*)1501 and DRB1(*)1502.

In Goodpasture syndrome, the yearly morbidity rate amounts to approx. 1 of 1.000.000 persons. The syndrome is characterised by bloody secretion and may lead to lung siderosis. Without treatment Goodpasture's syndrome has a very poor prognosis. Since the disease progresses rapidly, early diagnosis is imperative and can be confirmed by the detection of anti-GBM antibodies. A regular, quantitative determination of the autoantibody titer supports the monitoring of the therapy success, since the titer correlates with the disease activity.

Literature references

1. Bielsa I. **Update of systemic vasculitides nomenclature. International Chapel Hill Consensus Conference, 2012.** [Article in English, Spanish]. *Actas Dermosifiliogr* 106 (2015) 605-608.
2. Borza DB. **Autoepitopes and alloepitopes of type IV collagen: role in the molecular pathogenesis of anti-GBM antibody glomerulonephritis.** *Nephron Exp Nephrol* 106 (2007) 37-43.
3. EUROIMMUN AG. Stöcker W, Schlumberger W, Krüger C. **Alle Beiträge zum Thema Autoimmundiagnostik.** In: Gressner A, Arndt T. (Hrsg.). *Lexikon der Medizinischen Laboratoriumsdiagnostik*. 2. Auflage. Springer Medizin Verlag, Heidelberg (2012).
4. Hellmark T, Segelmark M. **Diagnosis and classification of Goodpasture's disease (anti-GBM).** *J Autoimmun* 48-49 (2014) 108-112.
5. Jennette JC, Falk RJ, Bacon PA, Basu N, Cid MC, Ferrario F, Flores-Suarez LF, Gross WL, Guillevin L, Hagen EC, Hoffman GS, Jayne DR, Kallenberg CG, Lamprecht P, Langford CA, Luqmani RA, Mahr AD, Matteson EL, Merkel PA, Ozen S, Pusey CD, Rasmussen N, Rees AJ, Scott DG, Specks U, Stone JH, Takahashi K, Watts RA. **2012 revised International Chapel Hill Consensus Conference Nomenclature of Vasculitides.** *Arthritis Rheum* 65 (2013) 1-11.

6. Sinico RA, Radice A, Corace C, Sabadini E, Bollini B. **Anti-glomerular basement membrane antibodies in the diagnosis of Goodpasture syndrome: a comparison of different assays.** Nephrol Dial Transplant 21 (2006) 397-401.
7. Yang R, Hellmark T, Zhao J, Cui Z, Segelmark M, Zhao MH, Wang HY. **Antigen and epitope specificity of anti-glomerular basement membrane antibodies in patients with goodpasture disease with or without anti-neutrophil cytoplasmic antibodies.** J Am Soc Nephrol 18 (2007) 1338-1343.
8. Silvariño R, Noboa O, Cervera R. **Anti-Glomerular Basement Membrane Antibodies.** IMAJ 16 (2014) 727-732.